

### **ORDERING INFORMATION**

2082 - Intelect® RPW unit - 100 V 2073 - Intelect® RPW unit - 120 V 2074 - Intelect® RPW unit - 230 V

#### **Standard Accessories:**

28700 - Standard accessory kit that includes the following:
28699 - D-ACTOR® hand piece applicator
28736 - Projectile (x3)
17638 - R15 15 mm ESWT transmitter
28724 - D20-S D-ACTOR® 20 mm transmitter
28737 - Sealing set
28738 - Guide tube (x2)
28739 - Cleaning brush

4248 - Conductor™ transmission gel 250 ml (8.5 oz) bottle (x1) 27465 - Patient data card (x1)

28668 - User manual on CD

### **Optional Accessories:**

28725 - D35 D-ACTOR® 35 mm transmitter 28726 - F15 Focus-Lens 15 mm transmitter 28728 - DI15 Deep Impact® 15 mm transmitter 28729 - C15 CERAma-x® ESWT 15 mm transmitter

28730 - V-ACTOR® hand piece applicator kit that includes the following:
28740 - V25 V-ACTOR® 25 mm vibration transmitter
28741 - V40 V-ACTOR® 40 mm vibration transmitter

DJO International
7 The Pines Business Park
Broad Street
Guildford
Surrey
GU3 3BH

Tel: +44 (0)1483 459 659 Fax: +44 (0)1483 459 470 www.DJOglobal.eu

**ENGLAND** 









Shockwave by Chattanooga $^{\mathsf{m}}$ 





# **INTELECT® RPW**

Intelect® RPW incorporates touch screen technology to ensure a high degree of simplicity. The easy-to-use menu-driven user interface guarantees reliable selection of all necessary parameters for treatment set-up as well as during patient treatment.

#### **USER PROFILE**

This Radial Pressure Wave clinical unit is ideally suited to the following applications:

- Orthopaedics
- Rheumatology
- Physical Medicine
- Physiotherapy
- Sport Clubs
- Acupuncture
- Aesthetic
- Veterinary

### **INDICATIONS**

- Myofascial Trigger Points (MTrP)
  - Localising and deactivating trigger points
- Activation of muscle and connective tissue
  - Increasing circulation
- Pulse vibration massage
- Disorder of tendon insertions
  - Plantar Fasciitis, heel pain, or heel spur
  - Tendinosis calcarea/supraspinatus-tendon
  - Radial and ulnar humeral epicondylitis
  - Achillodynia
  - Retropatellar pain syndrome
  - Tibial edge syndrome
  - Proximal iliotibial band friction syndrome/trochanteric
  - Insertional tendonitis



# **FEATURES & BENEFITS**

### **1) 2 CHANNEL OUTPUT CONNECTIONS**

Intelect® RPW allows the connection of 2 hand pieces for quick and efficient use of complementary applicator transmitters during patient treatment. 2 hand piece applicators available:

### **D-ACTOR®** Applicator

## Applications:

- Myofascial trigger points
- Other disorders of tendon insertions
- Activation of muscle and connective tissue
- Acupuncture shock wave therapy
- Cellulite treatment (stages I III)
- Muscle and connective tissue tightening
- Elimination of lymphatic congestion
- Body shaping
- Regulation of muscle tone
- Medical massage



#### V-ACTOR® Applicator

The V-ACTOR® applicator is an ideal complement to Radial Wave Therapy and contributes to long-lasting treatment results.

## Vibration treatment:

- Stimulates metabolism and accelerates the elimination of waste products
- Tightens tissue through muscular stimulation
- Improves skin firmness to restore a natural, younger appearance
- Offers a relaxing wellness experience

## Applications:

- Body shaping
- Lymphatic drainage
- Wellness massage
- Muscle and connective tissue activation



#### 2 STATE-OF-THE-ART NAVIGATION

This feature allows access to the desired functions at the touch of a finger. Treatment setup is quick and easy via Quick-link Indications, Clinical Protocols™, User Protocols or manual adjustment:

- 200 user-defined protocol storage slots
- Over 80 clinical protocols on board











# **3 UNIQUE OPTIMAL ENERGY LEVEL ADJUSTMENT**

- Unique ramp feature that allows gradual ramping of the intensity during treatment
- Provides an automatic and smooth intensity increase in tune with your patient
- Patient pain threshold can be recorded and treatment can be managed at the optimal painless energy level













## **4** CLINICAL RESOURCES

Intelect® RPW includes a full suite of reference resources available at the push of a button.

These resources include:

- Full colour graphics library of anatomical images and pathologies
- Hand piece placement images
- Radial pressure wave rationale
- Radial pressure wave contraindications



## **ADDITIONAL FEATURES & BENEFITS**

- Convenient storage platform allows secure and easy access of hand pieces, transmitters and gel
- 360° swivel rotation of interface
- Large storage drawer for storage of clinical supplies
- Multilingual software
- Patient Documentation of treatment sessions, pre and post treatment patient pain scales, pain mapping and pain types all captured on patient data cards







